

C1orf127 Antibody



PACO39942

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

C1orf127

Uniprot

Q8N9H9

Synonyms:

Uncharacterized protein C1orf127, C1orf127


Immunogen:

Recombinant Human Uncharacterized protein C1orf127 protein (401-648AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO39942 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic cancer using PACO39942 at dilution of 1:100.