

NANP Antibody, FITC conjugated



PACO39980

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytosol, N-acylneuraminate-9-phosphatase activity, cellular protein metabolic process, dephosphorylation, dolichol-linked oligosaccharide biosynthetic process, N-acetylneuraminate biosynthetic process, post-translational protein modification, protein N-linked glycosylation via asparagine.

Gene ID:

NANP

Uniprot

Q8TBE9

Synonyms:

N-acylneuraminate-9-phosphatase (EC 3.1.3.29) (Haloacid, dehalogenase-like hydrolase domain-containing protein 4) (Neu5Ac-9-Pase), NANP, C20orf147 HDHD4

Immunogen:

Recombinant Human N-acylneuraminate-9-phosphatase protein (1-248AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A