

# DUSP19 Antibody, HRP conjugated



PACO39991

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Has a dual specificity toward Ser/Thr and Tyr-containing proteins.

**Gene ID:**

DUSP19

**Uniprot**

Q8WTR2

**Synonyms:**

Dual specificity protein phosphatase 19 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity phosphatase TS-DSP1) (Low molecular weight dual specificity phosphatase 3) (LMW-DSP3) (Protein phosphatase SKRP1) (Stress-activated protein kinase pathway-regulating phosphatase 1) (SAPK pathway-regulating phosphatase 1), DUSP19, DUSP17, LMWDSP3, SKRP1

**Immunogen:**

Recombinant Human Dual specificity protein phosphatase 19 protein (1-217AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A