

CHMP7 Antibody, FITC conjugated



PACO40000

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays a role in the endosomal sorting pathway.

Gene ID:

CHMP7

Uniprot

Q8WUX9

Synonyms:

Charged multivesicular body protein 7 (Chromatin-modifying protein 7), CHMP7

Immunogen:

Recombinant Human Charged multivesicular body protein 7 protein (1-290AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A