

DCAKD Antibody, Biotin conjugated



PACO40009

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

membrane, dephospho-CoA kinase activity, coenzyme A biosynthetic process, phosphorylation.

Gene ID:

DCAKD

Uniprot

Q8WVC6

Synonyms:

Dephospho-CoA kinase domain-containing protein, DCAKD

Immunogen:

Recombinant Human Dephospho-CoA kinase domain-containing protein (1-231AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A