

TEKT4 Antibody



PACO40014

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Structural component of ciliary and flagellar microtubules. Forms filamentous polymers in the walls of ciliary and flagellar microtubules.

Gene ID:

TEKT4

Uniprot

Q8WW24

Synonyms:

Tektin-4, TEKT4

Immunogen:

Recombinant Human Tektin-4 protein (1-265AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A