

# MAIP1 Antibody, FITC conjugated



PACO40020

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MAIP1

**Uniprot**

Q8WWC4

**Synonyms:**

m-AAA protease-interacting protein 1, mitochondrial (Matrix AAA peptidase-interacting protein 1), MAIP1, C2orf47

**Immunogen:**

Recombinant Human m-AAA protease-interacting protein 1, mitochondrial protein (97-291AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A