

SLC7A3 Antibody



PACO40034

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Mediates the uptake of the cationic amino acid, arginine, lysine and ornithine in a sodium-independent manner.

Gene ID:

SLC7A3

Uniprot

Q8WY07

Synonyms:

Cationic amino acid, transporter 3 (CAT-3) (CAT3) (Cationic amino acid, transporter y+) (Solute carrier family 7 member 3), SLC7A3, ATRC3 CAT3


Immunogen:

Recombinant Human Cationic amino acid, transporter 3 protein (428-475AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of MCF-7 cells using PACO40034 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).