

# POMGNT1 Antibody, HRP conjugated



PACO40039

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Participates in O-mannosyl glycosylation. May be responsible for the synthesis of the GlcNAc(beta1-2)Man(alpha1-)O-Ser/Thr moiety on alpha-dystroglycan and other O-mannosylated proteins. Is specific for alpha linked terminal mannose and does not have MGAT3, MGAT4, MGAT5, MGAT7 or MGAT8 activity.

**Gene ID:**

POMGNT1

**Uniprot**

Q8WZA1

**Synonyms:**

Protein O-linked-mannose beta-1,2-N-acetylglucosaminyltransferase 1 (POMGnT1) (EC 2.4.1. -) (UDP-GlcNAc: alpha-D-mannoside beta-1,2-N-acetylglucosaminyltransferase I.2) (GnT I.2), POMGNT1, MGAT1.2

**Immunogen:**

Recombinant Human Protein O-linked-mannose beta-1,2-N-acetylglucosaminyltransferase 1 protein (59-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A