

FAM3C Antibody, HRP conjugated



PACO40047

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in retinal laminar formation. Promotes epithelial to mesenchymal transition.

Gene ID:

FAM3C

Uniprot

Q92520

Synonyms:

Protein FAM3C (Interleukin-like EMT inducer), FAM3C, ILEI

Immunogen:

Recombinant Human Protein FAM3C protein (31-227AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A