

# Pectate lyase 1 Antibody, HRP conjugated



PACO40063

---

## Product Information

**Size:**

50ug

**Reactivity:**

Chamaecyparis obtusa

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Has pectate lyase activity.

**Gene ID:****Uniprot**

Q96385

**Synonyms:**

Pectate lyase 1 (EC 4.2.2.2) (Major pollen allergen Cha o 1) (allergen Cha o 1)

**Immunogen:**

Recombinant Chamaecyparis obtusa Pectate lyase 1 protein (22-354AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A