

Pectate lyase 1 Antibody, Biotin conjugated



PACO40065

Product Information

Size:

50ug

Reactivity:

Chamaecyparis obtusa

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Has pectate lyase activity.

Gene ID:**Uniprot**

Q96385

Synonyms:

Pectate lyase 1 (EC 4.2.2.2) (Major pollen allergen Cha o 1) (allergen Cha o 1)

Immunogen:

Recombinant Chamaecyparis obtusa Pectate lyase 1 protein (22-354AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A