

TENT5B Antibody




PACO40070

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	TENT5B
Source:	Uniprot
Rabbit	Q96A09
Isotype:	Synonyms:
IgG	Putative nucleotidyltransferase FAM46B (EC 2.7.7. -), FAM46B
Applications:	Immunogen:
ELISA, WB	Recombinant Human Terminal nucleotidyltransferase 5B protein (2-425AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, WB:1:500-1:2000	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: TENT5B antibody at 6 μ g/ml. Lane 1: 293T whole cell lysate. Lane 2: HeLa whole cell lysate. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 47 kDa. Observed band size: 47 kDa.