

# PBXIP1 Antibody, HRP conjugated



PACO40083

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Regulator of pre-B-cell leukemia transcription factors (BPXs) function. Inhibits the binding of PBX1-HOX complex to DNA and blocks the transcriptional activity of E2A-PBX1. Tethers estrogen receptor-alpha (ESR1) to microtubules and allows them to influence estrogen receptors-alpha signaling.

**Gene ID:**

PBXIP1

**Uniprot**

Q96AQ6

**Synonyms:**

Pre-B-cell leukemia transcription factor-interacting protein 1 (Hematopoietic PBX-interacting protein), PBXIP1, HPIP

**Immunogen:**

Recombinant Human Pre-B-cell leukemia transcription factor-interacting protein 1 protein (361-600AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A