

PBXIP1 Antibody, FITC conjugated



PACO40084

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Regulator of pre-B-cell leukemia transcription factors (BPXs) function. Inhibits the binding of PBX1-HOX complex to DNA and blocks the transcriptional activity of E2A-PBX1. Tethers estrogen receptor-alpha (ESR1) to microtubules and allows them to influence estrogen receptors-alpha signaling.

Gene ID:

PBXIP1

Uniprot

Q96AQ6

Synonyms:

Pre-B-cell leukemia transcription factor-interacting protein 1 (Hematopoietic PBX-interacting protein), PBXIP1, HPIP

Immunogen:

Recombinant Human Pre-B-cell leukemia transcription factor-interacting protein 1 protein (361-600AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A