

C3orf49 Antibody



PACO40090

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

C3orf49

Uniprot

Q96BT1

Synonyms:

Putative uncharacterized protein C3orf49, C3orf49

Immunogen:

Recombinant Human Putative uncharacterized protein C3orf49 protein (1-292AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: C3orf49 antibody at 5 μ g/ml. Lane 1: U251 whole cell lysate. Lane 2: MCF-7 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 34 kDa. Observed band size: 34 kDa.