

GNPNAT1 Antibody, HRP conjugated



PACO40115

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytosol, glucosamine 6-phosphate N-acetyltransferase activity, identical protein binding, cellular protein metabolic process, dolichol-linked oligosaccharide biosynthetic process, post-translational protein modification, protein N-linked glycosylation via asparagine, UDP-N-acetylglucosamine biosynthetic process.

Gene ID:

GNPNAT1

Uniprot

Q96EK6

Synonyms:

Glucosamine 6-phosphate N-acetyltransferase (EC 2.3.1.4) (Phosphoglucosamine acetylase) (Phosphoglucosamine transacetylase), GNPNAT1, GNA1

Immunogen:

Recombinant Human Glucosamine 6-phosphate N-acetyltransferase protein (1-184AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A