

# GNPNAT1 Antibody, FITC conjugated



PACO40116

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

cytosol, glucosamine 6-phosphate N-acetyltransferase activity, identical protein binding, cellular protein metabolic process, dolichol-linked oligosaccharide biosynthetic process, post-translational protein modification, protein N-linked glycosylation via asparagine, UDP-N-acetylglucosamine biosynthetic process.

**Gene ID:**

GNPNAT1

**Uniprot**

Q96EK6

**Synonyms:**

Glucosamine 6-phosphate N-acetyltransferase (EC 2.3.1.4) (Phosphoglucosamine acetylase) (Phosphoglucosamine transacetylase), GNPNAT1, GNA1

**Immunogen:**

Recombinant Human Glucosamine 6-phosphate N-acetyltransferase protein (1-184AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A