

# PRAC1 Antibody



PACO40166

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PRAC1

**Uniprot**

Q96KF2

**Synonyms:**

Small nuclear protein PRAC1 (Prostate cancer susceptibility candidate protein 1) (Prostate, rectum and colon expressed gene protein), PRAC1, C17orf92 PRAC

**Immunogen:**

Recombinant Human Small nuclear protein PRAC1 protein (1-57AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A