

FAM118B Antibody



PACO40214

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

May play a role in Cajal bodies formation.

Gene ID:

FAM118B

Uniprot

Q9BPY3

Synonyms:

Protein FAM118B, FAM118B


Immunogen:

Recombinant Human Protein FAM118B protein (83-351AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: FAM118B antibody at 12 μ g/ml. Lane 1: Jurkat whole cell lysate. Lane 2: U251 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 39 kDa. Observed band size: 39 kDa.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO40214 at dilution of 1:100.