

S1PR3 Antibody, HRP conjugated



PACO40279

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Receptor for the lysosphingolipid sphingosine 1-phosphate (S1P). S1P is a bioactive lysophospholipid that elicits diverse physiological effect on most types of cells and tissues. When expressed in rat HTC4 hepatoma cells, is capable of mediating S1P-induced cell proliferation and suppression of apoptosis.

Gene ID:

S1PR3

Uniprot

Q99500

Synonyms:

Sphingosine 1-phosphate receptor 3 (S1P receptor 3) (S1P3) (Endothelial differentiation G-protein coupled receptor 3) (Sphingosine 1-phosphate receptor Edg-3) (S1P receptor Edg-3), S1PR3, EDG3

Immunogen:

Recombinant Human Sphingosine 1-phosphate receptor 3 protein (303-378AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A