

VP3 Antibody



PACO40290

Product Information

Size:

50ug

Reactivity:

Chicken anemia virus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

May act as transcriptional regulator. Induces apoptosis in infected cells. Element of infectious replication cycle.

Gene ID:

VP3

Uniprot

Q99152

Synonyms:

Apoptin, VP3


Immunogen:

Recombinant Chicken anemia virus Apoptin protein (1-121AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: Chicken anemia virus apoptin antibody at 2 μ g/ml + GST-Apoptin fusion protein at 0.4 μ g. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 40 kDa. Observed band size: 40 kDa.