

# VP3 Antibody, HRP conjugated



PACO40291

---

## Product Information

**Size:**

50ug

**Reactivity:**

Chicken anemia virus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May act as transcriptional regulator. Induces apoptosis in infected cells. Element of infectious replication cycle.

**Gene ID:**

VP3

**Uniprot**

Q99152

**Synonyms:**

Apoptin, VP3

**Immunogen:**

Recombinant Chicken anemia virus Apoptin protein (1-121AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A