

RIT2 Antibody, Biotin conjugated



PACO40313

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Binds and exchanges GTP and GDP.

Gene ID:

RIT2

Uniprot

Q99578

Synonyms:

GTP-binding protein Rit2 (Ras-like protein expressed in neurons) (Ras-like without CAAX protein 2), RIT2, RIN ROC2

Immunogen:

Recombinant Human GTP-binding protein Rit2 protein (1-150AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A