

RHOXF2 Antibody, Biotin conjugated



PACO40329

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

RHOXF2

Uniprot

Q9BQY4

Synonyms:

Rhox homeobox family member 2 (Paired-like homeobox protein PEPP-2) (Testis homeobox gene 1), RHOXF2, PEPP2 THG1

Immunogen:

Recombinant Human Rhox homeobox family member 2 protein (1-244AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A