

SDF4 Antibody, FITC conjugated



PACO40332

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May regulate calcium-dependent activities in the endoplasmic reticulum lumen or post-ER compartment.

Gene ID:

SDF4

Uniprot

Q9BRK5

Synonyms:

45 kDa calcium-binding protein (Cab45) (Stromal cell-derived factor 4) (SDF-4), SDF4, CAB45

Immunogen:

Recombinant Human 45 kDa calcium-binding protein (37-362AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A