

PACO40346

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Chaperone protein which promotes assembly of the 20S proteasome. May cooperate with PSMG1-PSMG2 heterodimers to orchestrate the correct assembly of proteasomes.

Gene ID:

PSMG3

Uniprot

Q9BT73

Synonyms:

Proteasome assembly chaperone 3 (PAC-3) (hPAC3), PSMG3, C7orf48 PAC3

Immunogen:

Recombinant Human Proteasome assembly chaperone 3 protein (1-122AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A