

EDEM2 Antibody, FITC conjugated



PACO40356

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Initiates the endoplasmic reticulum-associated degradation (ERAD) that targets misfolded glycoproteins for degradation in an N-glycan-dependent manner. Catalyzes the first mannose trimming step, from Man9GlcNAc2 to Man8GlcNAc2. Extracts misfolded glycoproteins, but not glycoproteins undergoing productive folding, from the calnexin cycle.

Gene ID:

EDEM2

Uniprot

Q9BV94

Synonyms:

ER degradation-enhancing alpha-mannosidase-like protein 2, EDEM2, C20orf31, C20orf49

Immunogen:

Recombinant Human ER degradation-enhancing alpha-mannosidase-like protein 2 protein (22-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A