

# STK33 Antibody, HRP conjugated



PACO40367

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Serine/threonine protein kinase which phosphorylates VIME. May play a specific role in the dynamic behavior of the intermediate filament cytoskeleton by phosphorylation of VIME. Not essential for the survival of KRAS-dependent AML cell lines.

**Gene ID:**

STK33

**Uniprot**

Q9BYT3

**Synonyms:**

Serine/threonine-protein kinase 33 (EC 2.7.11.1), STK33

**Immunogen:**

Recombinant Human Serine/threonine-protein kinase 33 protein (231-514AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A