

STK33 Antibody, FITC conjugated



PACO40368

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Serine/threonine protein kinase which phosphorylates VIME. May play a specific role in the dynamic behavior of the intermediate filament cytoskeleton by phosphorylation of VIME. Not essential for the survival of KRAS-dependent AML cell lines.

Gene ID:

STK33

Uniprot

Q9BYT3

Synonyms:

Serine/threonine-protein kinase 33 (EC 2.7.11.1), STK33

Immunogen:

Recombinant Human Serine/threonine-protein kinase 33 protein (231-514AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A