

ZDHHC5 Antibody, HRP conjugated



PACO40379

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Palmitoyl acyltransferase for the G-protein coupled receptor SSTR5. Also palmitoylates FLOT2.

Gene ID:

ZDHHC5

Uniprot

Q9C0B5

Synonyms:

Palmitoyltransferase ZDHHC5 (EC 2.3.1.225) (Zinc finger DHHC domain-containing protein 5) (DHHC-5) (Zinc finger protein 375), ZDHHC5, KIAA1748 ZNF375

Immunogen:

Recombinant Human Palmitoyltransferase ZDHHC5 protein (231-500AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A