

# ZDHHC5 Antibody, Biotin conjugated



PACO40381

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Palmitoyl acyltransferase for the G-protein coupled receptor SSTR5. Also palmitoylates FLOT2.

**Gene ID:**

ZDHHC5

**Uniprot**

Q9C0B5

**Synonyms:**

Palmitoyltransferase ZDHHC5 (EC 2.3.1.225) (Zinc finger DHHC domain-containing protein 5) (DHHC-5) (Zinc finger protein 375), ZDHHC5, KIAA1748 ZNF375

**Immunogen:**

Recombinant Human Palmitoyltransferase ZDHHC5 protein (231-500AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A