

GT-1 Antibody, FITC conjugated



PACO40384

Product Information

Size:

50ug

Reactivity:

Arabidopsis thaliana

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable transcription factor that binds specifically to the core DNA sequence 5'-GGTTAA-3'. May act as a molecular switch in response to light signals.

Gene ID:

GT-1

Uniprot

Q9FX53

Synonyms:

Trihelix transcription factor GT-1 (Trihelix DNA-binding protein GT-1), GT-1

Immunogen:

Recombinant Arabidopsis thaliana Trihelix transcription factor GT-1 protein (1-406AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A