

# SCAND2P Antibody



PACO40386

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

nucleus, sequence-specific DNA binding, transcription factor activity, sequence-specific DNA binding, regulation of transcription, DNA-templated.

**Gene ID:**

SCAND2P

**Uniprot**

Q9GZW5

**Synonyms:**

Putative SCAN domain-containing protein SCAND2P (SCAN domain-containing protein 2 pseudogene), SCAND2P, SCAND2

**Immunogen:**

Recombinant Human Putative SCAN domain-containing protein SCAND2P protein (1-152AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A