

# C11orf68 Antibody, HRP conjugated



PACO40399

---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	poly(A) RNA binding.
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	C11orf68
<b>Source:</b>	<b>Uniprot</b>
Rabbit	Q9H3H3
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	UPF0696 protein C11orf68 (Basophilic leukemia-expressed protein Bles03) (Protein p5326), C11orf68, BLES03
<b>Applications:</b>	<b>Immunogen:</b>
ELISA	Recombinant Human UPF0696 protein C11orf68 protein (1-251AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A