

# SERHL2 Antibody



PACO40406

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probable serine hydrolase. May be related to cell muscle hypertrophy.

**Gene ID:**

SERHL2

**Uniprot**

Q9H4I8

**Synonyms:**

Serine hydrolase-like protein 2 (EC 3.1. -. -), SERHL2, SERHL

**Immunogen:**

Recombinant Human Serine hydrolase-like protein 2 protein (1-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A