

ZDHHC6 Antibody, HRP conjugated



PACO40415

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Palmitoylates calnexin (CALX), which is required for its association with the ribosome-translocon complex and efficient folding of glycosylated proteins.

Gene ID:

ZDHHC6

Uniprot

Q9H6R6

Synonyms:

Palmitoyltransferase ZDHHC6 (EC 2.3.1.225) (Transmembrane protein H4) (Zinc finger DHHC domain-containing protein 6) (DHHC-6) (Zinc finger protein 376), ZDHHC6, ZNF376

Immunogen:

Recombinant Human Palmitoyltransferase ZDHHC6 protein (227-409AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A