

# ZDHHC6 Antibody, FITC conjugated



PACO40416

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Palmitoylates calnexin (CALX), which is required for its association with the ribosome-translocon complex and efficient folding of glycosylated proteins.

**Gene ID:**

ZDHHC6

**Uniprot**

Q9H6R6

**Synonyms:**

Palmitoyltransferase ZDHHC6 (EC 2.3.1.225) (Transmembrane protein H4) (Zinc finger DHHC domain-containing protein 6) (DHHC-6) (Zinc finger protein 376), ZDHHC6, ZNF376

**Immunogen:**

Recombinant Human Palmitoyltransferase ZDHHC6 protein (227-409AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A