

METTL7A Antibody, HRP conjugated



PACO40431

Product Information

Size:	Protein Background:
50ug	Probable methyltransferase.
Reactivity:	Gene ID:
Human	METTL7A
Source:	Uniprot
Rabbit	Q9H8H3
Isotype:	Synonyms:
IgG	Methyltransferase-like protein 7A (EC 2.1.1. -) (Protein AAM-B), METTL7A
Applications:	Immunogen:
ELISA	Recombinant Human Methyltransferase-like protein 7A protein (30-244AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A