

EFCAB1 Antibody



PACO40438

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

EFCAB1

Uniprot

Q9HAE3

Synonyms:

EF-hand calcium-binding domain-containing protein 1, EFCAB1

Immunogen:

Recombinant Human EF-hand calcium-binding domain-containing protein 1 protein (1-211AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A