

HSPA6 Antibody




PACO40466

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	HSPA6
Source:	Uniprot
Rabbit	P17066
Isotype:	Synonyms:
IgG	Heat shock 70 kDa protein 6 (Heat shock 70 kDa protein B'), HSPA6, HSP70B'
Applications:	Immunogen:
ELISA, WB	Recombinant Human Heat shock 70 kDa protein 6 protein (1-637AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, WB:1:1000-1:5000	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: HSPA6 antibody at 2 μ g/ml. Lane 1: MCF-7 whole cell lysate. Lane 2: A431 whole cell lysate. Lane 3: Jurkat whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 72 kDa. Observed band size: 72 kDa.