

DBNDD2 Antibody



PACO40470

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May modulate the activity of casein kinase-1. Inhibits CSNK1D autophosphorylation.

Gene ID:

DBNDD2

Uniprot

Q9BQY9

Synonyms:

Dysbindin domain-containing protein 2 (Casein kinase-1 binding protein) (CK1BP) (HSMNP1), DBNDD2, C20orf35


Immunogen:

Recombinant Human Dysbindin domain-containing protein 2 protein (1-161AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO40470 at dilution of 1:100.