

ADAT1 Antibody, Biotin conjugated



PACO40489

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Specifically deaminates adenosine-37 to inosine in tRNA-Ala.

Gene ID:

ADAT1

Uniprot

Q9BUB4

Synonyms:

tRNA-specific adenosine deaminase 1 (hADAT1) (EC 3.5.4.34) (tRNA-specific adenosine-37 deaminase), ADAT1

Immunogen:

Recombinant Human tRNA-specific adenosine deaminase 1 protein (1-501AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A