

SNRNP25 Antibody, HRP conjugated



PACO40495

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytoplasm, intercellular bridge, nucleoplasm, nucleus, U12-type spliceosomal complex, gene expression, mRNA splicing, via spliceosome, RNA splicing.

Gene ID:

SNRNP25

Uniprot

Q9BV90

Synonyms:

U11/U12 small nuclear ribonucleoprotein 25 kDa protein (U11/U12 snRNP 25 kDa protein) (U11/U12-25K) (Minus-99 protein), SNRNP25, C16orf33

Immunogen:

Recombinant Human U11/U12 small nuclear ribonucleoprotein 25 kDa protein (1-132AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A