

# SNRNP25 Antibody, FITC conjugated



PACO40496

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

cytoplasm, intercellular bridge, nucleoplasm, nucleus, U12-type spliceosomal complex, gene expression, mRNA splicing, via spliceosome, RNA splicing.

**Gene ID:**

SNRNP25

**Uniprot**

Q9BV90

**Synonyms:**

U11/U12 small nuclear ribonucleoprotein 25 kDa protein (U11/U12 snRNP 25 kDa protein) (U11/U12-25K) (Minus-99 protein), SNRNP25, C16orf33

**Immunogen:**

Recombinant Human U11/U12 small nuclear ribonucleoprotein 25 kDa protein (1-132AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A