

TTLL2 Antibody



PACO40502

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Probable tubulin polyglutamylase that forms polyglutamate side chains on tubulin. Probably acts when complexed with other proteins.

Gene ID:

TTLL2

Uniprot

Q9BWW7

Synonyms:

Probable tubulin polyglutamylase TTLL2 (EC 6) (Testis-specific protein NYD-TSPG) (Tubulin--tyrosine ligase-like protein 2), TTLL2, C6orf104


Immunogen:

Recombinant Human Probable tubulin polyglutamylase TTLL2 protein (401-592AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO40502 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic cancer using PACO40502 at dilution of 1:100.