

KRTAP3-1 Antibody



PACO40514

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

In the hair cortex, hair keratin intermediate filaments are embedded in an interfilamentous matrix, consisting of hair keratin-associated proteins (KRTAP), which are essential for the formation of a rigid and resistant hair shaft through their extensive disulfide bond cross-linking with abundant cysteine residues of hair keratins. The matrix proteins include the high-sulfur and high-glycine-tyrosine keratins.

Gene ID:

KRTAP3-1

Uniprot

Q9BYR8

Synonyms:

Keratin-associated protein 3-1 (High sulfur keratin-associated protein 3.1) (Keratin-associated protein 3.1), KRTAP3-1, KAP3.1 KRTAP3.1

Immunogen:

Recombinant Human Keratin-associated protein 3-1 protein (2-98AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A