

# FXVD6 Antibody, HRP conjugated



PACO40539

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

integral component of plasma membrane, plasma membrane, sodium channel regulator activity, ion transmembrane transport, regulation of sodium ion transmembrane transporter activity, transmembrane transport.

**Gene ID:**

FXVD6

**Uniprot**

Q9H0Q3

**Synonyms:**

FXVD domain-containing ion transport regulator 6 (Phosphohippolin), FXVD6

**Immunogen:**

Recombinant Human FXVD domain-containing ion transport regulator 6 protein (59-95AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A