

FXVD6 Antibody, FITC conjugated



PACO40540

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

integral component of plasma membrane, plasma membrane, sodium channel regulator activity, ion transmembrane transport, regulation of sodium ion transmembrane transporter activity, transmembrane transport.

Gene ID:

FXVD6

Uniprot

Q9H0Q3

Synonyms:

FXVD domain-containing ion transport regulator 6 (Phosphohippolin), FXVD6

Immunogen:

Recombinant Human FXVD domain-containing ion transport regulator 6 protein (59-95AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A