

LMAN2L Antibody, FITC conjugated



PACO40544

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in the regulation of export from the endoplasmic reticulum of a subset of glycoproteins. May function as a regulator of ERGIC-53.

Gene ID:

LMAN2L

Uniprot

Q9H0V9

Synonyms:

VIP36-like protein (Lectin mannose-binding 2-like) (LMAN2-like protein), LMAN2L, VIPL

Immunogen:

Recombinant Human VIP36-like protein (45-313AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A